



# Objectives-Lenses



## Objective (Lenses), Soachale



For standard eye-pieces, with fixing-screw.

Simple handling by spring-lock for take-off the eye-piece.

With c-mount-thread to the camera head.

MM-25-Z-028	Focusing Lens F = 28 mm
MM-25-Z-035	Focusing Lens F = 35 mm
MM-25-Z-040	Focusing Lens F = 40 mm
MM-25-Z-100	Vario Focusing Lens F = 25 40 mm
MM-25-Z-001	Adaptor-plastics, connection between c-mount and 17 mm thread

